



LIVEY 950BT ANC WIRELESS BLUETOOTH HEADSET

Wireless Headset for On-the-Go and
Business Professionals





FOR OFFICE & ENTERPRISES, ON-THE-GO

Crystal clear calls with dual mic ENC & Advanced Noise Cancellation (ANC) with Intelligent Audio and seamless communication with low-latency.



Triple Mode
Connectivity Options



Hybrid Active Noise
Cancellation (ANC)



Dual-ENC Noise
Cancelling Mic



MS Teams Optimized BT
Dongle (USB-C) + Wireless

Compatible





**BT for stable and
reliable connectivity**



**(ANC) technology for
noise reduction**



**Comfortable
Ergonomic design**



**High-Quality
Leather Cushions**

KEY FEATURES

- Dual-mic Environmental Noise Cancellation (ENC) for clear calls
- Triple Mode Connectivity Options
- MS Teams Optimized BT Dongle (USB-C) + Wireless
- Advanced Bluetooth Technology, Bluetooth 5.4 for stable and reliable connectivity
- Hybrid Active Noise Cancellation (ANC) technology for noise reduction
- 2.4Ghz Wireless frequency when used with Type-C dongle
- Also can be used as wired mode using the 3.5mm Stereo Cable
- Long-Lasting Battery Life, Up to 100 hours of playback time (ANC OFF) and 60 hours (ANC ON) at 50% volume
- High-Quality Audio sensitivity
- Comes with Comfortable Ergonomic design

TECHNICAL SPECIFICATION

Pairing name	LIVEY 950BT Headset
BT Version	BT 5.4
Support	A2DP/AVRCP/HSP/HFP
Operation Range	10m
Rated power	10mW (one ear)
Max power	30mW (one ear)
Frequency range	20-20KHz
Impedance	36Ω±15%
Driver Sensitivity	110dB±3dB
Distortion Rate	≤4.0% (1kHz/IEC318)
Driver size	Φ40mm
Dimensions	197*188*79 mm
Weight	290g

Battery	Lithium-ion Polymer Battery
Capacity	1000mAh
Charging Time	less than 3hrs
ANC	Hybrid, -38dB
ANC Mic Sensitivity	-31dB, SNR:66dB
ANC Mic Impedance	2.2kΩ
ENC	Dual-mic ENC Call Noise Cancellation
ENC Mic Sensitivity	-31dB, SNR:68dB
ENC Mic Impedance	2.2kΩ
Latency (Dongle Mode)	less than 30ms
Buttons	2 (ANC button + MFB)
Indicator Light	Red and blue dual-color LED working indicator light, ANC green light
Play time (ANC OFF)	100 hrs (50% volume)
Play time (ANC ON)	60 hrs (50% volume)